



Hain Lifescience – a Bruker Company

## Tuberculosis

### ● Rapid diagnostics of tuberculosis and its resistances

Our molecular test systems allow for an efficient TB screening and subsequent detection of MDR- or XDR-TB. **FluoroType® MTB** is the innovative test system for every lab and enables the direct detection of *M. tuberculosis* (MTB) complex. A resistance testing is possible with **GenoType MTBDRplus** and **GenoType MTBDRsl** – the globally established tests for MDR- and XDR-TB.

**TBCheck MPT64** is a robust test for the detection of *M. tuberculosis* complex directly from culture. After confirmation of MTB complex a PCR-based differentiation is possible with **GenoType MTBC**.

### Your benefits with the TB product series from Hain Lifescience

- Fast and reliable results
- Step-wise diagnostics
- Cost efficiency
- CE-marked

## TB Product Series

TB is the most prevalent infectious disease worldwide. A single patient with active TB may infect 10 to 15 other people each year. Four parameters are of crucial importance for the control of the disease:

- Early diagnostics
- Prevention of the spread of the disease
- Effective treatment with antitubercotics
- Prevention of the development of drug resistances

The TB product series from Hain Lifescience offers rapid, easy and cost-efficient diagnostic systems that are the prerequisite for an effective treatment and confinement of tuberculosis.

### TB screening

Early and reliable diagnostics are the basis for a specific and thus successful tuberculosis treatment. Culture methods are time-consuming and laborious. In contrast, nucleic acid amplification tests, which allow for a fast screening have proven themselves in practice. The **FluoroType® MTB** test uses PCR and an innovative fluorescent-based technology for the detection of *M. tuberculosis* complex directly from pulmonary and extrapulmonary clinical specimens. The results are available within only three hours –therefore **FluoroType® MTB** is the ideal TB screening test.

### Drug susceptibility testing (DST)

The increase of Multidrug-resistant (MDR-)TB is an alarming and ongoing global issue. MDR-TB is defined as TB that is resistant to at least rifampicin and isoniazid, the two most powerful first-line drugs. In order to prevent the further spread of resistant TB and offer the most appropriate therapy rapid and direct detection of MDR-TB is mandatory. **GenoType MTBDRplus** is based on PCR and the DNA•STRIP technology and allows the detection of *M. tuberculosis* complex and its resistance against rifampicin and isoniazid directly from clinical specimens. Extensively drug-resistant (XDR-)TB is defined as MDR-TB with further resistance to fluoroquinolones and a second-line agent (amikacin, kanamycin or capreomycin). Diagnostics and treatment of XDR-TB are even more challenging as those strains leave patients nearly without any treatment options. **GenoType MTBDRsl** can be performed subsequently to **GenoType MTBDRplus** using the same DNA isolate – thus an efficient testing for XDR-TB is possible.

### Culture Identification


Culture confirmation still plays a considerable role in TB diagnostics. The **TBCheck MPT64** assay detects and confirms *M. tuberculosis* complex via the MPT64 antigen from culture material. Depending on the results further diagnostics can be initiated efficiently.

### Differentiation

In former times the routinely performed differentiation of MTB complex was only possible with phenotypic and biochemical methods which are laborious and time-consuming. **GenoType MTBC** is a PCR test for the differentiation of MTB complex from cultivated specimens.

## FluoroType® MTB

Fluorescence-based molecular genetic assay for direct detection of *M. tuberculosis* complex



### Characteristics of FluoroType® MTB

The **FluoroType® MTB** permits the molecular genetic identification of the *M. tuberculosis* complex directly from pulmonary and extrapulmonary patient specimens via the detection of the IS6110 element. The assay is rapid, simple and highly sensitive.

### Test principle of FluoroType® MTB

**FluoroType® MTB** is based on PCR and the **FluoroType®** technology. Mycobacterial DNA is extracted from the patient specimen and specifically amplified via PCR. Then fluorescence-labelled probes are bound to single-stranded amplicons. Changes in fluorescence intensity are measured and displayed as a melting curve. The evaluation is done by the test-specific software. Amplification and detection run fully automated in the **FluoroCycler®** instruments. This ensures maximum user-friendliness and efficient diagnostics with reliable results at one glance.


### Benefits of using FluoroType® MTB

- **Fast:** The test system provides reliable results within three hours. This allows for an important time advantage in TB diagnostics.
- **User-friendly:** A ready-to-use amplification mix already containing the Taq polymerase is provided with the kit. Amplification and detection run fully automated in a closed system. Evaluation and result interpretation are done by the test-specific software. The test procedure is very simple and perfectly suitable for low to mid throughput.
- **Flexible:** DNA extraction can be performed either manually with **FluoroLyse** or automated using the nucleic acid extraction instrument **GenoXtract®**.



# GenoType MTBDR<sub>plus</sub> VER 2.0

Molecular genetic assay for detection of *M. tuberculosis* complex and its resistances to rifampicin and/or isoniazid



## Characteristics of GenoType MTBDR<sub>plus</sub>

The **GenoType MTBDR<sub>plus</sub>** enables the simultaneous molecular genetic identification of

- the *M. tuberculosis* complex
- its resistance to rifampicin by the detection of the most common mutations in the *rpoB* gene
- its resistance to isoniazid (For the detection of high level isoniazid resistance the *katG* gene and for low level isoniazid resistance the promoter region of the *inhA* gene is examined.)

from smear-positive or -negative pulmonary clinical specimens or cultivated samples.

## Test principle of GenoType MTBDR<sub>plus</sub>

**GenoType MTBDR<sub>plus</sub>** is based on PCR and the DNA•STRIP technology. Mycobacterial DNA is extracted from the patient specimen or cultivated material, specifically amplified via PCR and detected on a membrane strip using reverse hybridization and an enzymatic colour reaction. Valid results are documented by internal controls, Conjugate and Amplification Control.


## Benefits of using GenoType MTBDR<sub>plus</sub>

- **Efficient:** *M. tuberculosis* complex and its resistances to rifampicin and isoniazid are simultaneously detected in a single patient specimen. The test is therefore perfectly suitable for MDR-TB screening, for the identification of MTB complex and mono-resistances. Pulmonary patient specimens and cultivated samples can be used as starting material.
- **Rapid:** Results are available within five hours compared to several months with conventional DST.
- **User-friendly:** A ready-to-use amplification mix including the Taq polymerase is provided with the kit.
- **Flexible:** DNA extraction can be performed either manually or automated using the nucleic acid isolation instrument GenoXtract®. Amplification, detection and evaluation can also be automated. The test is thus suitable for low, mid and high throughput.
- **Cost-efficient:** For the implementation only minimum technical equipment is required, therefore an economical set-up is possible for laboratories of every potential size.



# GenoType MTBDR<sub>s/</sub>

Molecular genetic assay for detection of *M. tuberculosis* complex and its resistances to fluoroquinolones and aminoglycosides/cyclic peptides (and ethambutol)



	GenoType MTBDR <sub>s/</sub> VER 1.0	GenoType MTBDR <sub>s/</sub> VER 2.0
Detection of	<i>M. tuberculosis</i> complex and its resistances to fluoroquinolones, aminoglycosides/cyclic peptides and ethambutol	<i>M. tuberculosis</i> complex and its resistances to fluoroquinolones and aminoglycosides/cyclic peptides
Sample Material	smear-positive pulmonary and cultivated samples	smear-positive and -negative pulmonary and cultivated samples
Ethambutol	Mutations in the <i>embB</i> gene that are involved in ethambutol resistance	
	✓	-
Fluoroquinolone	Mutations in the <i>gyrB</i> gene that are involved in fluoroquinolone resistance	
	-	✓
Kanamycin	Mutations in the <i>eis</i> gene that are involved in kanamycin resistance	
	-	✓

## Test principle of GenoType MTBDR<sub>s/</sub>


**GenoType MTBDR<sub>s/</sub>** is based on PCR and the DNA•STRIP technology. Mycobacterial DNA is extracted from the patient specimen or cultivated material, specifically amplified via PCR and detected on a membrane strip using reverse hybridization and an enzymatic colour reaction.

## Benefits of using GenoType MTBDR<sub>s/</sub>

- **Sensitive detection:** The first version of **GenoType MTBDR<sub>s/</sub>** can be processed from smear-positive pulmonary or cultivated samples. The second version is even more sensitive and can therefore also be performed using smear-negative pulmonary samples.
- **Efficient diagnosis:** Both test systems are perfectly suitable for the detection of XDR-TB in patients previously diagnosed with MDR-TB. For step-wise diagnostics the test systems can be performed subsequent to **GenoType MTBDR<sub>plus</sub>** using the same DNA isolate.
- **Rapid results:** Results are available within five hours in comparison to several weeks when using conventional methods.

# TBCheck MPT64 VER 1.0

Rapid immunochromatographic assay for detection of *Mycobacterium tuberculosis* complex



## Characteristics of TBCheck MPT64

The **TBCheck MPT64** assay allows the identification of the MPT64 antigen from cultivated liquid samples. This antigen is highly specific for *M. tuberculosis* complex and thus suitable for its identification directly from culture. As the MPT64 antigen is only present in MTB complex subsequent discrimination from nontuberculous mycobacteria is also possible.

## Test principle of TBCheck MPT64


**TBCheck MPT64** is based on an immunochromatographic assay principle. A droplet of the positive culture is placed on the lateral flow strip. On the strip the secreted MPT64 antigens are marked with gold and migrate to a specific binding site. This reaction leads to an accumulation of gold at the binding site and subsequently to a visible band on the strip. The control area shows the efficiency of the gold binding – therefore, valid results are always guaranteed.

## Benefits of using TBCheck MPT64

- **Rapid detection:** **TBCheck MPT64** allows the rapid detection of *M. tuberculosis* complex and discrimination from NTM within 10 minutes. Therefore, rapid results are guaranteed and further testing is promptly possible.
- **Indication for further diagnostics:** The results of **TBCheck MPT64** enable a sound choice for further diagnostics. Depending on the result, differentiation of *M. tuberculosis* complex or NTM is indicated.
- **Confirmation:** The assay can be used to confirm *M. tuberculosis* complex before drug susceptibility testing is performed.

# GenoType MTBC VER 1.X

Molecular genetic assay for differentiation of *M. tuberculosis* complex



## Characteristics of GenoType MTBC

The **GenoType MTBC** permits the molecular genetic identification of *M. tuberculosis*, *M. bovis ssp. bovis*, *M. bovis ssp. caprae*, *M. africanum*, *M. microti*, *M. canettii* and the vaccine strain *M. bovis* BCG (Bacille Calmette-Guérin) from cultivated samples.

## Test principle of GenoType MTBC

**GenoType MTBC** is based on PCR and the DNA•STRIP technology. Mycobacterial DNA is extracted from cultivated material, specifically amplified via PCR and detected on a membrane strip using reverse hybridization and an enzymatic colour reaction. Valid results are documented by the Conjugate Control. The Universal Control displays the presence of mycobacteria and gram-positive bacteria with high G+C content. The MTBC control shows that members of the MTB complex are present.

## Benefits of using GenoType MTBC

- **Efficient:** Simultaneous detection and differentiation of species belonging to the *M. tuberculosis* complex with a single processing. As starting material solid or liquid cultivated material can be used.
- **Rapid:** The results are available within five hours compared to several weeks with conventional methods.
- **Reliable:** Internal controls document valid results and thus ensure high diagnostic reliability.

## Mycobacteria Product Series

### TB screening

#### FluoroType® MTB VER 1.0

Detection of *M. tuberculosis* complex from patient specimens

### Drug susceptibility testing

#### FluoroType® MTBDR VER 2.0

Single-tube detection of *M. tuberculosis* complex and its resistances to rifampicin and isoniazid from patient specimens or cultures

#### GenoType MTBDRplus VER 2.0

Detection of *M. tuberculosis* complex and its resistances to rifampicin and isoniazid from patient specimens or cultures

#### GenoType MTBDRsl VER 1.0

Detection of *M. tuberculosis* complex and its resistances to fluoroquinolones, aminoglycosides/cyclic peptides and ethambutol from patient specimens or cultures

#### GenoType MTBDRsl VER 2.0

Detection of *M. tuberculosis* complex and its resistances to fluoroquinolones and aminoglycoside/cyclic peptides from patient specimens or cultures

### Differentiation

#### GenoType MTBC VER 1.X

Differentiation of *M. tuberculosis* complex from cultures

#### GenoType CMdirect VER 1.0

Detection of *M. tuberculosis* complex and more than 20 clinically relevant NTM from patient specimens

#### GenoType Mycobacterium CM VER 2.0

Detection of *M. tuberculosis* complex and more than 20 clinically relevant NTM from cultures

#### GenoType Mycobacterium AS VER 1.0

Detection of 19 further NTM from cultures

### Differentiation and drug susceptibility testing

#### GenoType NTM-DR VER 1.0

Detection of important NTM and their resistances to aminoglycosides and macrolides from cultures

### Culture identification

#### TBCheck MPT64 VER 1.0

Rapid detection of *M. tuberculosis* complex from liquid cultures

### Leprosy

#### GenoType LepraeDR VER 1.0

Detection of *M. leprae* and its resistances to rifampicin, ofloxacin and dapsone from patient specimens

Please contact your local representative for availability in your country. Not for sale in the USA.



Hardwiesenstrasse 1  
72147 Nehren - Germany  
Phone +49 (0) 7473-9451-0



[info.mdx.de@bruker.com](mailto:info.mdx.de@bruker.com) - [www.hain-lifescience.de](http://www.hain-lifescience.de)