

Overview PCR programs for the test systems of Hain Lifescience

GenoType MRSA Direct	MRSA DIR
GenoType MRSA	MRSA HOT
GenoType Staphylococcus GenoType BC grampositive GenoType BC gramnegative micro-IDent® micro-IDent®<i>plus</i>	30CYCLES
GenoType CDiff	30CYCLES for culture samples; 40CYCLES for stool samples
GenoType PAI-1 GenoType HH GenoType AAT GenoType LCT	HOT 30; HOT 35 for GENO•CARD
ThromboType® ThromboType®<i>plus</i>	HOT 30; HOT 30 or HOT 35 for GENO•CARD
GenoType MTHFR	HOT 30; HOT 30 for GENO•CARD
GenoType ApoE	HOT 35; HOT 40 for GENO•CARD
GenoType HelicoDR	HOT 30 for culture samples; HOT 35 for biopsy samples
GenoType Mycobacteria Direct	MYCO DIR
GenoType MTBDR<i>plus</i> GenoType MTBDR<i>sl</i>	HOT 30 for culture samples; HOT 40 for direct samples
GenoType LepraeDR	HOT 40
GenoType MTBC GenoType Mycobacterium CM GenoType Mycobacterium AS GenoType EHEC GenoType Enterococcus	HOT 30
GenoQuick® MRSA	GQFST1
GenoQuick® HLA-B27	GQHOT 30; GQHOT 35 for GENO•CARD
GenoQuick® Bordetella	GQFST1
GenoQuick® CT	GQHOT 40
GenoQuick® MTB	GQMTB