


For a correctly sized printout, adjust Adobe Acrobat Reader's Print Menu settings as shown here:




Use these settings, print only page 2 and/or 3 of this document.

GenoType MTBDRsl 12


VER 1.0
00317-0315-02-3

dd	mm	yyyy

#			
---	--	---	---



TUB	FLQ	AG/CP	EMB
<i>gyrA</i> WT	sensitive	sensitive	sensitive
<i>gyrA</i> MUT	resistant	resistant	resistant
<i>rrs</i> WT	sensitive	sensitive	sensitive
<i>rrs</i> MUT	resistant	resistant	resistant
<i>embB</i> WT	sensitive	sensitive	sensitive
<i>embB</i> MUT	resistant	resistant	resistant



LOT _____
  HYB _____ min
 STR _____ min
 SUB _____ min
 _____


GenoType MTBDRsl 12

VER 2.0

00317A-0215-03-1

dd	mm	yyyy

# 		
--	---	---



TUB	<i>gyrA</i> WT	<i>gyrA</i> MUT	<i>gyrB</i> WT	<i>gyrB</i> MUT	<i>rrs</i> WT	<i>rrs</i> MUT	<i>eis</i> WT	<i>eis</i> MUT	Phenotypic resistance				
									FLQ	KAN/AMK/CAP	KAN/CAP/NO	KAN/AMK/CAP/NO	low-level KAN

#	Person	Person with envelope	TUB	<i>gyrA</i> WT	<i>gyrA</i> MUT	<i>gyrB</i> WT	<i>gyrB</i> MUT	<i>rrs</i> WT	<i>rrs</i> MUT	<i>eis</i> WT	<i>eis</i> MUT	FLQ	KAN/AMK/CAP	KAN/CAP/NO	KAN/AMK/CAP/NO	low-level KAN

LOT _____

 HYB _____ min

 STR _____ min

 SUB _____ min

 _____

